

Anti-CD49d antibody [HP2/1] (STJ16101920)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC/IP
Host / Source	Mouse
Reactivity	Human/Rat/Monkey/Bovine/Pig/Goat/Rabbit/Horse/Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---