

## Anti-CD79a antibody [HM57] (STJ16101926)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rabbit/Rat/Pig/Bovine/Guinea Pig/Monkey

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

### TARGET INFORMATION

Gene ID	<a href="#">973</a>
Gene Symbol	<a href="#">CD79A</a>
UniProt ID	<a href="#">CD79A_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---