

Anti-Chymotrypsin antibody [CHYMO10-104.5] (STJ16101930)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/ELISA/IF/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Formulation	PBS with <0.1% sodium azide
Isotype	IgG3
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---