

Anti-Insulin/Proinsulin antibody [D3E7] (STJ16101934)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-F/ELISA/IHC-P
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Bovine/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---