

## Anti-Lysozyme antibody [BGN/0696/5B1] (STJ16101935)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IHC-P
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Formulation	PBS with <0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---