

Anti-MAdCAM-1 antibody [17F5] (STJ16101936)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/FC/IHC-P/WB |
| Host / Source | Mouse |
| Reactivity | Human/Camel |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Purified |
| Formulation | PBS with <0.1% sodium azide |
| Isotype | IgG1 |
| Storage Instruction | Store aliquots at -20°C, and thawed aliquots at 2-8°C. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 8174 |
| Gene Symbol | MADCAM1 |
| UniProt ID | MADCA_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|