

Anti-N-Cadherin antibody [13A9] (STJ16101940)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IF/IP/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

TARGET INFORMATION

Gene ID	1000
Gene Symbol	CDH2
UniProt ID	CADH2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---