

## Anti-Parathyroid Hormone antibody [BGN/1F8] (STJ16101942)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/ELISA/IHC-P
Host / Source	Mouse
Reactivity	Human/Rat/Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

### TARGET INFORMATION

Gene ID	<a href="#">5741</a>
Gene Symbol	<a href="#">PTH</a>
UniProt ID	<a href="#">PTHY_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---