

## Anti-Pepsinogen I antibody [8003] (STJ16101943)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/IHC-P        |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 0.5 mg/mL  |
| Conjugation         | Unconjugated   |
| Purification        | Purified   |
| Formulation         | PBS with <0.1% sodium azide                            |
| Isotype             | IgG1   |
| Storage Instruction | Store aliquots at -20°C, and thawed aliquots at 2-8°C. |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|