

Anti-Synaptophysin antibody [SP15] (STJ16101945)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Hamster/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Formulation	PBS with <0.1% sodium azide
Isotype	IgM
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

TARGET INFORMATION

Gene ID	6855
Gene Symbol	SYP
UniProt ID	SYPH_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---