

Anti-TALIN-1 antibody [97H6] (STJ16101946)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/IP/IHC-P/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

TARGET INFORMATION

Gene ID	7094
Gene Symbol	TLN1
UniProt ID	TLN1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---