

Anti-Lassa Virus NP antibody [EBS-I-315] (STJ16101968)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/WB/IF/IHC/LF
Host / Source	Mouse
Reactivity	Lassa Virus NP

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: 0, 1-100 Mu g/ml; tracer: 0, 001-100 Mu g/ml for 30 min at RT). Immunoblotting (1-2 Mu g/ml). Immunofluorescence (0, 5-1, 0 Mu g/ml). Immunohistology (1-2 Mu g/ml for 30 min at RT). Lateral flow (solid phase: 0, 4 Mu g per cm; tra
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---