

Anti-MUC2 antibody [EBS-T-233] (STJ16101978)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/WB/IF/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow Cytometry (0, 5-1, 0 µg/million cells in 0, 1ml). Immunoblotting (1-2 µg/ml). Immunofluorescence (1-2 µg/ml). Immunohistology (formalin-fixed: 1-2 µg/ml for 30 min at RT; requires boiling tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 1
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4583
Gene Symbol	MUC2
UniProt ID	MUC2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---