

Anti-ASGPR/CD301 antibody [121A5.01] {Alexa-fluor-488} (STJ16101992)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FACS |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|----------------------------------|
| Clonality | Monoclonal |
| Concentration | 0.5 mg/mL |
| Conjugation | Alexa-fluor-488 |
| Purification | Affinity Purified |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG1 |
| Storage Instruction | Store at -20°C for up to 1 year. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|