

**Anti-Bet V1 antibody [H.G26F6B2.06] (STJ16102041)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/RIPA
Host / Source	Human
Reactivity	Pollen

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG3
Storage Instruction	Store at -20°C for up to 1 year.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---