

## Anti-Bet V1 antibody [H.G3A3C2.04] (STJ16102045)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/RIPA
Host / Source	Human
Reactivity	Pollen

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2
Storage Instruction	Store at -20°C for up to 1 year.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------