

## Anti-CD1a antibody [214A9.01] {Alexa-fluor-647} (STJ16102049)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year.

### TARGET INFORMATION

Gene ID	<a href="#">909</a>
Gene Symbol	<a href="#">CD1A</a>
UniProt ID	<a href="#">CD1A_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---