

Anti-CD14 antibody [MOP9.25] (STJ16102058)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FACS/Cell Depletion/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2b, k
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	929
Gene Symbol	CD14
UniProt ID	CD14_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---