

**Anti-CD40L/CD154 antibody [LL2] (STJ16102066)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC/IF/ELISA        |
| Host / Source | Mouse              |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                  |
|---------------------|----------------------------------|
| Clonality           | Monoclonal                       |
| Concentration       | 0.5 mg/mL                        |
| Conjugation         | Unconjugated                     |
| Purification        | Affinity Purified                |
| Formulation         | TRIS, PBS, pH7 and antibiotics   |
| Isotype             | IgG1                             |
| Storage Instruction | Store at -20°C for up to 1 year. |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|