

Anti-CD38 antibody [HB7] (STJ16102076)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FACS
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	952
Gene Symbol	CD38
UniProt ID	CD38_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---