

**Anti-DC-LAMP/CD208 antibody [104G4] {Biotin} (STJ16102089)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FACS/IHC-F
Host / Source	Mouse
Reactivity	Human/Sheep

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---