

**Anti-DC-LAMP/CD208 antibody [1010E1.01] {Alexa-fluor-546} (STJ16102092)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IF
Host / Source	Rat
Reactivity	Mouse/Dog/Cat/Human/Sheep/Mouse

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-546
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage Instruction	Store at-20°C for up to 1 year.

**ADDITIONAL INFORMATION**

Note	<b>STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.</b>
------	--