

Anti-DC-LAMP/CD208 antibody [1010E1.01] {Biotin} (STJ16102094)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FACS/IHC-F
Host / Source	Rat
Reactivity	Mouse/Dog/Cat/Human/Sheep/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---