

Anti-DC-LAMP/CD208 antibody [208B5] {Alexa-fluor-647} (STJ16102108)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FACS/IHC-F
Host / Source	Mouse
Reactivity	Mouse/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG3
Storage Instruction	Store at-20°C for up to 1 year.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---