

Anti-DC-SIGN/CD209 antibody [103G2.07] {Biotin} (STJ16102125)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FACS/IHC-F/IHC-P
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1, k
Storage Instruction	Store at -20°C for up to 1 year.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---