

**Anti-DC-SIGN/CD209 antibody [120C8] {Alexa-fluor-647} (STJ16102142)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FACS
Host / Source	Mouse
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2b, k
Storage Instruction	Store at -20°C for up to 1 year.

**ADDITIONAL INFORMATION**

Note	<b>STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.</b>
------	--