

Anti-DERP-1 antibody [H.AK2D] {HRPO} (STJ16102157)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Human
Reactivity	House Dust Mite

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	HRPO
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1k
Storage Instruction	Store at -20°C for up to 1 year.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---