

**Anti-DORA/IgSF6 antibody [104A10.01] {Alexa-fluor-488} (STJ16102159)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FACS/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Dog

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at-20°C for up to 1 year.

**ADDITIONAL INFORMATION**

Note	<b>STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.</b>
------	--