

Anti-DORA/IgSF6 antibody [104A10.01] {Alexa-fluor-647} (STJ16102161)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FACS/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at-20°C for up to 1 year.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--