

Anti-IL3R Alpha/CD123 antibody [106C2.02] {Alexa-fluor-647} (STJ16102193)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---