

## Anti-IL10 antibody [12G8] {Alexa-fluor-546} (STJ16102198)

### GENERAL INFORMATION

Product Type	Primary antibodies
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-546
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year.

### TARGET INFORMATION

Gene ID	<a href="#">16153</a>
Gene Symbol	<a href="#">IL10</a>
UniProt ID	<a href="#">IL10_MOUSE</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---