

**Anti-Langerin/CD207 antibody [929F3.01] {Alexa-fluor-488} (STJ16102246)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FACS
Host / Source	Rat
Reactivity	Human/Mouse/Rat/Pig

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---