

Anti-mouse Langerin/CD207 antibody [923B7] {Alexa-fluor-488} (STJ16102266)**GENERAL INFORMATION**

| | |
|---------------|---------------------|
| Product Type | Primary antibodies |
| Applications | IF/FACS |
| Host / Source | Rat |
| Reactivity | Human/Rat/Mouse/Pig |

PRODUCT PROPERTIES

| | |
|---------------------|---------------------------------|
| Clonality | Monoclonal |
| Concentration | 0.5 mg/mL |
| Conjugation | Alexa-fluor-488 |
| Purification | Affinity Purified |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG2a |
| Storage Instruction | Store at-20°C for up to 1 year. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|