

## Anti-mouse Langerin/CD207 antibody [923B7] {Alexa-fluor-546} (STJ16102267)

### GENERAL INFORMATION

|               |                     |
|---------------|---------------------|
| Product Type  | Primary antibodies  |
| Applications  | IF                  |
| Host / Source | Rat                 |
| Reactivity    | Human/Rat/Mouse/Pig |

### PRODUCT PROPERTIES

|                     |                                 |
|---------------------|---------------------------------|
| Clonality           | Monoclonal                      |
| Concentration       | 0.5 mg/mL                       |
| Conjugation         | Alexa-fluor-546                 |
| Purification        | Affinity Purified               |
| Formulation         | TRIS, PBS, pH7 and antibiotics  |
| Isotype             | IgG2a                           |
| Storage Instruction | Store at-20°C for up to 1 year. |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|