

Anti-mouse Langerin/CD207 antibody [923B7] {Alexa-fluor-647} (STJ16102268)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FACS
Host / Source	Rat
Reactivity	Human/Rat/Mouse/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage Instruction	Store at-20°C for up to 1 year.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---