

**Anti-MIP-3 Alpha/CCL20 antibody [206D9.05] (STJ16102303)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-P/IHC-F/ELISA
Host / Source	Mouse
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---