

Anti-Neuropilin/CD304 antibody [211H6.01] {Biotin} (STJ16102311)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IP/IHC-F/FACS
Host / Source	Mouse
Reactivity	Human/Dog/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---