

Anti-TLR3/CD283 antibody [619F7.06] {Biotin} (STJ16102326)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/FACS/WB
Host / Source	Mouse
Reactivity	Human/Dog/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	7098
Gene Symbol	TLR3
UniProt ID	TLR3_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------