

## Anti-TLR3/CD283 antibody [722E2.01] {Biotin} (STJ16102331)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB
Host / Source	Mouse
Reactivity	Human/Rat/Mouse/Dog/Sheep/Pig

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year.

### TARGET INFORMATION

Gene ID	<a href="#">7098</a>
Gene Symbol	<a href="#">TLR3</a>
UniProt ID	<a href="#">TLR3_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---