

**Anti-TLR3/CD283 antibody [1213F10] {Alexa-fluor-647} (STJ16102345)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/WB/IP/IF
Host / Source	Mouse
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

**TARGET INFORMATION**

Gene ID	<a href="#">7098</a>
Gene Symbol	<a href="#">TLR3</a>
UniProt ID	<a href="#">TLR3_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---