

**Anti-TLR8/CD288 antibody [303F1.14] {Alexa-fluor-488} (STJ16102355)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IF/FACS
Host / Source	Mouse
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year.

**TARGET INFORMATION**

Gene ID	<a href="#">170744</a>
Gene Symbol	<a href="#">Tlr8</a>
UniProt ID	<a href="#">TLR8_MOUSE</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---