

**Anti-TLR8/CD288 antibody [303F1.14] {Alexa-fluor-546} (STJ16102356)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IF                 |
| Host / Source | Mouse              |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                  |
|---------------------|----------------------------------|
| Clonality           | Monoclonal                       |
| Concentration       | 0.5 mg/mL                        |
| Conjugation         | Alexa-fluor-546                  |
| Purification        | Affinity Purified                |
| Formulation         | TRIS, PBS, pH7 and antibiotics   |
| Isotype             | IgG2a                            |
| Storage Instruction | Store at -20°C for up to 1 year. |

**TARGET INFORMATION**

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">170744</a>     |
| Gene Symbol | <a href="#">Tlr8</a>       |
| UniProt ID  | <a href="#">TLR8_MOUSE</a> |

**ADDITIONAL INFORMATION**

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|