

Anti-HIV-1 gp120/160 antibody [109A2.01] {Alexa-fluor-647} (STJ16102418)**GENERAL INFORMATION**

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FACS/IF |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---------------------------------|
| Clonality | Monoclonal |
| Concentration | 0.5 mg/mL |
| Conjugation | Alexa-fluor-647 |
| Purification | Affinity Purified |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG1 |
| Storage Instruction | Store at-20°C for up to 1 year. |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|