

**Anti-HIV-1 gp41 antibody [417D3.08/clone 2] {Alexa-fluor-647} (STJ16102427)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FACS/IF/Neut
Host / Source	Mouse
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at-20°C for up to 1 year.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---