

**Anti-HIV-1 gp41 antibody [425D11.12/clone 3] {Alexa-fluor-647} (STJ16102433)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/FACS/IF/Neut |
| Host / Source | Mouse              |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                 |
|---------------------|---------------------------------|
| Clonality           | Monoclonal                      |
| Concentration       | 0.5 mg/mL                       |
| Conjugation         | Alexa-fluor-647                 |
| Purification        | Affinity Purified               |
| Formulation         | TRIS, PBS, pH7 and antibiotics  |
| Isotype             | IgG1                            |
| Storage Instruction | Store at-20°C for up to 1 year. |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|