

Anti-Myd88 antibody [603E10.05] {Alexa-fluor-546} (STJ16102437)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FACS
Host / Source	Mouse
Reactivity	Human/Dog/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-546
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2b
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	17874
Gene Symbol	Myd88
UniProt ID	MYD88_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---