

## Anti-Myd88 antibody [603E10.05] (STJ16102439)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FACS
Host / Source	Mouse
Reactivity	Human/Dog/Pig

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2b
Storage Instruction	Store at -20°C for up to 1 year.

### TARGET INFORMATION

Gene ID	<a href="#">17874</a>
Gene Symbol	<a href="#">Myd88</a>
UniProt ID	<a href="#">MYD88_MOUSE</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---