

**Anti-Alpha 4 Beta 7 antibody [204C5] {Alexa-fluor-647} (STJ16102455)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FACS/IF
Host / Source	Mouse
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at-20°C for up to 1 year.

**ADDITIONAL INFORMATION**

Note	<b>STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.</b>
------	--