

## Anti-Human IgE antibody {Alexa-fluor-647} (STJ16102520)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF
Host / Source	Human
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgE
Storage Instruction	Store at -20°C for up to 1 year.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---