

## Anti-Rat IgG2a antibody {Alexa-fluor-647} (STJ16102545)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FACS/IF            |
| Host / Source | Rat                |
| Reactivity    | Rat                |

### PRODUCT PROPERTIES

|                     |                                  |
|---------------------|----------------------------------|
| Clonality           | Polyclonal                       |
| Concentration       | 0.5 mg/mL                        |
| Conjugation         | Alexa-fluor-647                  |
| Purification        | Affinity Purified                |
| Formulation         | TRIS, PBS, pH7 and antibiotics   |
| Isotype             | IgG2a                            |
| Storage Instruction | Store at -20°C for up to 1 year. |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|